

C1QTNF3 Antibody, HRP conjugated



PACO25345

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

C1QTNF3 is a secretory protein produced by chondrogenic precursors & proliferating chondrocytes, and belongs to a novel C1q family of proteins.

Gene ID:

C1QTNF3

Uniprot

Q9BXJ4

Synonyms:

Complement C1q tumor necrosis factor-related protein 3 (Collagenous repeat-containing sequence 26 kDa protein) (CORS26) (Secretory protein CORS26), C1QTNF3, CTRP3

Immunogen:

Recombinant Human Complement C1q tumor necrosis factor-related protein 3 protein (23-246AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A