

# C1QTNF3 Antibody, FITC conjugated



PACO25346

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

C1QTNF3 is a secretory protein produced by chondrogenic precursors & proliferating chondrocytes, and belongs to a novel C1q family of proteins.

**Gene ID:**

C1QTNF3

**Uniprot**

Q9BXJ4

**Synonyms:**

Complement C1q tumor necrosis factor-related protein 3 (Collagenous repeat-containing sequence 26 kDa protein) (CORS26) (Secretory protein CORS26), C1QTNF3, CTRP3

**Immunogen:**

Recombinant Human Complement C1q tumor necrosis factor-related protein 3 protein (23-246AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A