

EIF4A2 Antibody, FITC conjugated



PACO25366

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

EIF4A2

Uniprot

Q14240

Synonyms:

Eukaryotic initiation factor 4A-II (eIF-4A-II) (eIF4A-II) (EC 3.6.4.13) (ATP-dependent RNA helicase eIF4A-2), EIF4A2, DDX2B EIF4F

Immunogen:

Recombinant Human Eukaryotic initiation factor 4A-II protein (1-407AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A