

# EIF4A2 Antibody, Biotin conjugated



PACO25367

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

EIF4A2

**Uniprot**

Q14240

**Synonyms:**

Eukaryotic initiation factor 4A-II (eIF-4A-II) (eIF4A-II) (EC 3.6.4.13) (ATP-dependent RNA helicase eIF4A-2), EIF4A2, DDX2B EIF4F

**Immunogen:**

Recombinant Human Eukaryotic initiation factor 4A-II protein (1-407AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A