

COX5B Antibody, Biotin conjugated



PACO25399

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

This protein is one of the nuclear-coded polypeptide chains of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport.

Gene ID:

COX5B

Uniprot

P10606

Synonyms:

Cytochrome c oxidase subunit 5B, mitochondrial (Cytochrome c oxidase polypeptide Vb), COX5B

Immunogen:

Recombinant Human Cytochrome c oxidase subunit 5B, mitochondrial protein (32-129AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A