

ERCC6L2 Antibody, FITC conjugated



PACO25418

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in early DNA damage response.

Gene ID:

ERCC6L2

Uniprot

Q5T890

Synonyms:

DNA excision repair protein ERCC-6-like 2 (EC 3.6.4. -) (DNA repair and recombination protein RAD26-like), ERCC6L2, C9orf102 RAD26L

Immunogen:

Recombinant Human DNA excision repair protein ERCC-6-like 2 protein (72-371AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A