

# CA1 Antibody, FITC conjugated



PACO25428

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Reversible hydration of carbon dioxide. Can hydrates cyanamide to urea.

**Gene ID:**

CA1

**Uniprot**

P00915

**Synonyms:**

Carbonic anhydrase 1 (EC 4.2.1.1) (Carbonate dehydratase I) (Carbonic anhydrase B) (CAB) (Carbonic anhydrase I) (CA-I), CA1

**Immunogen:**

Recombinant Human Carbonic anhydrase 1 protein (2-261AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A