

# Ca1 Antibody, FITC conjugated



PACO25429

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Reversible hydration of carbon dioxide.

**Gene ID:**

Ca1

**Uniprot**

B0BNN3

**Synonyms:**

Carbonic anhydrase 1 (EC 4.2.1.1) (Carbonate dehydratase I) (Carbonic anhydrase I) (CA-I), Ca1, Car1

**Immunogen:**

Recombinant Rat Carbonic anhydrase 1 protein (2-261AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A