

Ca1 Antibody, FITC conjugated



PACO25434

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Reversible hydration of carbon dioxide.

Gene ID:

Ca1

Uniprot

B0BNN3

Synonyms:

Carbonic anhydrase 1 (EC 4.2.1.1) (Carbonate dehydratase I) (Carbonic anhydrase I) (CA-I), Ca1, Car1

Immunogen:

Recombinant Rat Carbonic anhydrase 1 protein (2-261AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A