

CA12 Antibody, HRP conjugated



PACO25441

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Reversible hydration of carbon dioxide.

Gene ID:

CA12

Uniprot

O43570

Synonyms:

Carbonic anhydrase 12 (EC 4.2.1.1) (Carbonate dehydratase XII) (Carbonic anhydrase XII) (CA-XII) (Tumor antigen HOM-RCC-3.1.3), CA12

Immunogen:

Recombinant Human Carbonic anhydrase 12 protein (25-301AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A