

# CA12 Antibody, HRP conjugated



PACO25445

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Reversible hydration of carbon dioxide.

**Gene ID:**

CA12

**Uniprot**

O43570

**Synonyms:**

Carbonic anhydrase 12 (EC 4.2.1.1) (Carbonate dehydratase XII) (Carbonic anhydrase XII) (CA-XII) (Tumor antigen HOM-RCC-3.1.3), CA12

**Immunogen:**

Recombinant Human Carbonic anhydrase 12 protein (25-301AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A