

PACO25464

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Potent and specific inhibitor of CaM-kinase II (CAMK2).

**Gene ID:**

Camk2n1

**Uniprot**

Q6QWF9

**Synonyms:**

Calcium/calmodulin-dependent protein kinase II inhibitor 1 (calcium/calmodulin-dependent protein kinase II inhibitor alpha) (mCaMKIINalpha), Camk2n1

**Immunogen:**

Recombinant Mouse Calcium/calmodulin-dependent protein kinase II inhibitor 1 protein (1-78AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A