

Camk2n1 Antibody, FITC conjugated



PACO25466

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Potent and specific inhibitor of CaM-kinase II (CAMK2).

Gene ID:

Camk2n1

Uniprot

Q6QWF9

Synonyms:

Calcium/calmodulin-dependent protein kinase II inhibitor 1 (calcium/calmodulin-dependent protein kinase II inhibitor alpha) (mCaMKIINalpha), Camk2n1

Immunogen:

Recombinant Mouse Calcium/calmodulin-dependent protein kinase II inhibitor 1 protein (1-78AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A