

# CAPN2 Antibody, FITC conjugated



PACO25480

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Calcium-regulated non-lysosomal thiol-protease which catalyze limited proteolysis of substrates involved in cytoskeletal remodeling and signal transduction.

**Gene ID:**

CAPN2

**Uniprot**

P17655

**Synonyms:**

Calpain-2 catalytic subunit (EC 3.4.22.53) (Calcium-activated neutral proteinase 2) (CANP 2) (Calpain M-type) (Calpain large polypeptide L2) (Calpain-2 large subunit) (Millimolar-calpain) (M-calpain), CAPN2, CANPL2

**Immunogen:**

Recombinant Human Calpain-2 catalytic subunit protein (20-700AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A