

CAPN2 Antibody, Biotin conjugated



PACO25483

Product Information

Size:

50ug

Reactivity:

Pig

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Calcium-regulated non-lysosomal thiol-protease which catalyze limited proteolysis of substrates involved in cytoskeletal remodeling and signal transduction. Proteolytically cleaves MYOC at 'Arg-226'.

Gene ID:

CAPN2

Uniprot

P43367

Synonyms:

Calpain-2 catalytic subunit (EC 3.4.22.53) (Calcium-activated neutral proteinase 2) (CANP 2) (Calpain M-type) (Calpain-2 large subunit) (Millimolar-calpain) (M-calpain) (Fragment), CAPN2

Immunogen:

Recombinant Pig Calpain-2 catalytic subunit protein (1-324AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A