

CARS Antibody, HRP conjugated



PACO25493

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytoplasm, cytosol, ATP binding, cysteine-tRNA ligase activity, protein homodimerization activity, tRNA binding, cysteinyl-tRNA aminoacylation, gene expression, tRNA aminoacylation for protein translation.

Gene ID:

CARS

Uniprot

P49589

Synonyms:

Cysteine--tRNA ligase, cytoplasmic (EC 6.1.1.16) (Cysteinyl-tRNA synthetase) (CysRS), CARS

Immunogen:

Recombinant Human Cysteine--tRNA ligase, cytoplasmic protein (201-500AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A