

# HSBP1 Antibody



PACO25508

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Negative regulator of the heat shock response. Negatively affects HSF1 DNA-binding activity. May have a role in the suppression of the activation of the stress response during the aging process.

**Gene ID:**

HSBP1

**Uniprot**

O75506

**Synonyms:**

Heat shock factor-binding protein 1 (Nasopharyngeal carcinoma-associated antigen 13) (NPC-A-13), HSBP1, HSF1BP

**Immunogen:**


Recombinant Human Heat shock factor-binding protein 1 protein (1-75AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human lung tissue using PACO25508 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO25508 at dilution of 1:100.