

Hemoglobin Antibody, Biotin conjugated



PACO25555

Product Information

Size:

50ug

Reactivity:

Guinea Pig

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

Uniprot

Synonyms:

Hemoglobin

Immunogen:

Guinea Pig Hemoglobin protein (Native Protein).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A