

CD160 Antibody, FITC conjugated



PACO25582

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Receptor showing broad specificity for both classical and non-classical MHC class I molecules.

Gene ID:

CD160

Uniprot

O95971

Synonyms:

CD160 antigen (Natural killer cell receptor BY55) (CD antigen CD160), CD160, BY55

Immunogen:

Recombinant Human CD160 antigen protein (27-159AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A