

Cd82 Antibody, FITC conjugated



PACO25638

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Associates with CD4 or CD8 and delivers costimulatory signals for the TCR/CD3 pathway.

Gene ID:

Cd82

Uniprot

P40237

Synonyms:

CD82 antigen (C33 antigen) (IA4) (Inducible membrane protein R2) (Metastasis suppressor Kangai-1 homolog) (CD antigen CD82), Cd82, Kai1

Immunogen:

Recombinant Mouse CD82 antigen protein (111-227AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A