

# CDKN2AIPNL Antibody



PACO25688

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

CDKN2AIPNL

**Uniprot**

Q96HQ2

**Synonyms:**

CDKN2AIP N-terminal-like protein (CDKN2A-interacting protein N-terminal-like protein), CDKN2AIPNL

**Immunogen:**


Recombinant Human CDKN2AIP N-terminal-like protein (1-116AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded mouse liver tissue using PACO25688 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO25688 at dilution of 1:100.