

CEACAM4 Antibody, Biotin conjugated



PACO25708

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Granulocyte orphan receptor that acts as a trigger efficient phagocytosis of attached particles.

Gene ID:

CEACAM4

Uniprot

O75871

Synonyms:

Carcinoembryonic antigen-related cell adhesion molecule 4 (Carcinoembryonic antigen CGM7) (Non-specific cross-reacting antigen W236), CEACAM4, CGM7

Immunogen:

Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 4 protein (36-155AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A