

CEACAM7 Antibody, FITC conjugated



PACO25711

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

integral component of membrane, plasma membrane.

Gene ID:

CEACAM7

Uniprot

Q14002

Synonyms:

Carcinoembryonic antigen-related cell adhesion molecule 7 (Carcinoembryonic antigen CGM2), CEACAM7, CGM2

Immunogen:

Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 7 protein (36-242AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A