

Cela2a Antibody, HRP conjugated



PACO25714

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts upon elastin.

Gene ID:

Cela2a

Uniprot

P05208

Synonyms:

Chymotrypsin-like elastase family member 2A (EC 3.4.21.71) (Elastase-2) (Elastase-2A), Cela2a, Ela-2 Ela2 Ela2a

Immunogen:

Recombinant Mouse Chymotrypsin-like elastase family member 2A protein (31-271AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A