

CHRNA1 Antibody, Biotin conjugated



PACO25764

Product Information

Size:

50ug

Reactivity:

Torpedo californica

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

After binding acetylcholine, the AChR responds by an extensive change in conformation that affects all subunits and leads to opening of an ion-conducting channel across the plasma membrane.

Gene ID:

CHRNA1

Uniprot

P02710

Synonyms:

Acetylcholine receptor subunit alpha, CHRNA1

Immunogen:

Recombinant Torpedo californica Acetylcholine receptor subunit α protein (25-234AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A