

CNBP Antibody, Biotin conjugated



PACO25832

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Single-stranded DNA-binding protein, with specificity to the sterol regulatory element (SRE). Involved in sterol-mediated repression.

Gene ID:

CNBP

Uniprot

P62633

Synonyms:

Cellular nucleic acid, binding protein (CNBP) (Zinc finger protein 9), CNBP, RNF163 ZNF9

Immunogen:

Recombinant Human Cellular nucleic acid, binding protein (2-170AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A