

CNKSR2 Antibody, HRP conjugated



PACO25834

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May function as an adapter protein or regulator of Ras signaling pathways.

Gene ID:

CNKSR2

Uniprot

Q8WXI2

Synonyms:

Connector enhancer of kinase suppressor of ras 2 (Connector enhancer of KSR 2) (CNK homolog protein 2) (CNK2), CNKSR2, CNK2 KIAA0902 KSR2

Immunogen:

Recombinant Human Connector enhancer of kinase suppressor of ras 2 protein (650-800AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A