

COL12A1 Antibody, HRP conjugated



PACO25854

Product Information

Size:

50ug

Reactivity:

Chicken

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Type XII collagen interacts with type I collagen-containing fibrils, the COL1 domain could be associated with the surface of the fibrils, and the COL2 and NC3 domains may be localized in the periferibrillar matrix.

Gene ID:

COL12A1

Uniprot

P13944

Synonyms:

Collagen alpha-1(XII) chain (Fibrochimerin), COL12A1

Immunogen:

Recombinant Chicken Collagen alpha-1(XII) chain protein (139-311AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A