

Col1a1 Antibody



PACO25857

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:10-1:100

Protein Background:

Gene ID:

Col1a1

Uniprot

P02454

Synonyms:

Collagen alpha-1(I) chain (Alpha-1 type I collagen), Col1a1


Immunogen:

Recombinant Rat Collagen alpha-1(I) chain protein (955-1207AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO25857 at dilution of 1:20.