

# Col1a1 Antibody, Biotin conjugated



PACO25860

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	
<b>Reactivity:</b>	<b>Gene ID:</b>
Rat	Col1a1
<b>Source:</b>	<b>Uniprot</b>
Rabbit	P02454
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Collagen alpha-1(I) chain (Alpha-1 type I collagen), Col1a1
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Rat Collagen alpha-1(I) chain protein (955-1207AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A