

COPB1 Antibody, Biotin conjugated



PACO25876

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

COPB1

Uniprot

P53618

Synonyms:

Coatomer subunit beta (Beta-coat protein) (Beta-COP), COPB1, COPB

Immunogen:

Recombinant Human Coatomer subunit β protein (2-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A