

EMC8 Antibody, HRP conjugated



PACO25886

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytoplasm, ER membrane protein complex, membrane, mitochondrion, nucleus.

Gene ID:

EMC8

Uniprot

O43402

Synonyms:

ER membrane protein complex subunit 8 (Neighbor of COX4) (Protein FAM158B), EMC8, C16orf2 C16orf4 COX4AL COX4NB FAM158B NOC4

Immunogen:

Recombinant Human ER membrane protein complex subunit 8 protein (1-210AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A