

CPA6 Antibody, FITC conjugated



PACO25899

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in the proteolytic inactivation of enkephalins and neurotensin in some brain areas. May convert inactive angiotensin I into the biologically active angiotensin II.

Gene ID:

CPA6

Uniprot

Q8N4T0

Synonyms:

Carboxypeptidase A6 (EC 3.4.17. -), CPA6, CPAH

Immunogen:

Recombinant Human Carboxypeptidase A6 protein (130-437AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A