

RPL8 Antibody, Biotin conjugated



PACO25928

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

This protein can be partially incorporated into E.coli polysomes in vivo, indicating it can replace the endogenous protein.

Gene ID:

RPL8

Uniprot

P62917

Synonyms:

60S ribosomal protein L8 (Large ribosomal subunit protein uL2), RPL8

Immunogen:

Recombinant Human 60S ribosomal protein L8 protein (3-257AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A