

# CRHR1 Antibody, HRP conjugated



PACO25930

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CRHR1

**Uniprot**

P34998

**Synonyms:**

Corticotropin-releasing factor receptor 1 (CRF-R-1) (CRF-R1) (CRFR-1) (Corticotropin-releasing hormone receptor 1) (CRH-R-1) (CRH-R1), CRHR1, CRFR CRFR1 CRHR

**Immunogen:**

Recombinant Human Corticotropin-releasing factor receptor 1 protein (24-121AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A