

# CRIP2 Antibody, HRP conjugated



PACO25938

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CRIP2

**Uniprot**

P52943

**Synonyms:**

Cysteine-rich protein 2 (CRP-2) (Protein ESP1), CRIP2, CRP2

**Immunogen:**

Recombinant Human Cysteine-rich protein 2 protein (1-208AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A