

CRIP2 Antibody, FITC conjugated



PACO25939

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	CRIP2
Source:	Uniprot
Rabbit	P52943
Isotype:	Synonyms:
IgG	Cysteine-rich protein 2 (CRP-2) (Protein ESP1), CRIP2, CRP2
Applications:	Immunogen:
ELISA	Recombinant Human Cysteine-rich protein 2 protein (1-208AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A