

CRIP2 Antibody, Biotin conjugated



PACO25944

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in the cytoskeletal anchoring of DLG4 in excitatory synapses.

Gene ID:

CRIP2

Uniprot

Q9P021

Synonyms:

Cysteine-rich PDZ-binding protein (Cysteine-rich interactor of PDZ three) (Cysteine-rich interactor of PDZ3), CRIP2

Immunogen:

Recombinant Human Cysteine-rich PDZ-binding protein (1-101AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A