

# CRYBA4 Antibody, FITC conjugated



PACO25959

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Crystallins are the dominant structural components of the vertebrate eye lens.

**Gene ID:**

CRYBA4

**Uniprot**

P53673

**Synonyms:**

Beta-crystallin A4 (Beta-A4 crystallin), CRYBA4

**Immunogen:**

Recombinant Human Beta-crystallin A4 protein (2-196AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A