

CRYZL1 Antibody, Biotin conjugated



PACO25964

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

cytosol, NADP binding, NADPH: quinone reductase activity, quinone metabolic process.

Gene ID:

CRYZL1

Uniprot

O95825

Synonyms:

Quinone oxidoreductase-like protein 1 (EC 1. -. -. -) (Protein 4P11) (Quinone oxidoreductase homolog 1) (QOH-1) (Zeta-crystallin homolog), CRYZL1, 4P11

Immunogen:

Recombinant Human Quinone oxidoreductase-like protein 1 protein (1-87AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A