

Csad Antibody, FITC conjugated



PACO25971

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

carboxy-lyase activity, L-cysteine catabolic process to taurine, taurine metabolic process.

Gene ID:

Csad

Uniprot

Q9DBE0

Synonyms:

Cysteine sulfinic acid, decarboxylase (EC 4.1.1.29) (Cysteine-sulfinic acid decarboxylase) (Sulfinic acid decarboxylase), Csad

Immunogen:

Recombinant Mouse Cysteine sulfinic acid, decarboxylase protein (1-493AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A