

# CSF1R Antibody



PACO25977

---

## Product Information

**Size:**

50ug

**Reactivity:**

Felis silvestris catus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CSF1R

**Uniprot**

P13369

**Synonyms:**

Macrophage colony-stimulating factor 1 receptor (CSF-1 receptor) (CSF-1-R) (CSF-1R) (M-CSF-R) (EC 2.7.10.1) (Proto-oncogene c-Fms) (CD antigen CD115), CSF1R, FMS

**Immunogen:**

Recombinant Felis silvestris catus Macrophage colony-stimulating factor 1 receptor protein (20-514AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A