

CSNK1D Antibody, FITC conjugated



PACO25999

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

CSNK1D

Uniprot

P48730

Synonyms:

Casein kinase I isoform delta (CKI-delta) (CKId) (EC 2.7.11.1) (Tau-protein kinase CSNK1D) (EC 2.7.11.26), CSNK1D, HCKID

Immunogen:

Recombinant Human Casein kinase I isoform delta protein (152-398AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A