

Ctsk Antibody, HRP conjugated



PACO26034

Product Information

Size:

50ug

Reactivity:

Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Closely involved in osteoclastic bone resorption and may participate partially in the disorder of bone remodeling. Displays potent endoprotease activity against fibrinogen at acid, pH. May play an important role in extracellular matrix degradation By similarity.

Gene ID:

Ctsk

Uniprot

O35186

Synonyms:

Cathepsin K (EC 3.4.22.38), Ctsk

Immunogen:

Recombinant Rat Cathepsin K protein (115-329AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A