

CXorf1 Antibody



PACO26053

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

CXorf1

Uniprot

O96002

Synonyms:

Putative transmembrane protein CXorf1, CXorf1


Immunogen:

Recombinant Human Putative transmembrane protein CXorf1 protein (79-111AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO26053 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human brain tissue using PACO26053 at dilution of 1:100.