

# HYPM Antibody



PACO26061

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

nuclear chromatin, DNA binding, chromatin organization, chromatin silencing.

**Gene ID:**

HYPM

**Uniprot**

O75409

**Synonyms:**

Huntingtin-interacting protein M (Huntingtin yeast partner M), HYPM, CXorf27

**Immunogen:**


Recombinant Human Huntingtin-interacting protein M protein (1-117AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human colon cancer using PACO26061 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human pancreatic cancer using PACO26061 at dilution of 1:100.