

# CYB5B Antibody, FITC conjugated



PACO26071

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Cytochrome b5 is a membrane bound hemoprotein which function as an electron carrier for several membrane bound oxygenases.

**Gene ID:**

CYB5B

**Uniprot**

O43169

**Synonyms:**

Cytochrome b5 type B (Cytochrome b5 outer mitochondrial membrane isoform), CYB5B, CYB5M OMB5

**Immunogen:**

Recombinant Human Cytochrome b5 type B protein (2-118AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A