

CYB5B Antibody, Biotin conjugated



PACO26072

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cytochrome b5 is a membrane bound hemoprotein which function as an electron carrier for several membrane bound oxygenases.

Gene ID:

CYB5B

Uniprot

O43169

Synonyms:

Cytochrome b5 type B (Cytochrome b5 outer mitochondrial membrane isoform), CYB5B, CYB5M OMB5

Immunogen:

Recombinant Human Cytochrome b5 type B protein (2-118AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A