

CYGB Antibody, Biotin conjugated



PACO26088

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May have a protective function during conditions of oxidative stress. May be involved in intracellular oxygen storage or transfer.

Gene ID:

CYGB

Uniprot

Q8WWM9

Synonyms:

Cytoglobin (Histoglobin) (HGb) (Stellate cell activation-associated protein), CYGB, STAP

Immunogen:

Recombinant Human Cytoglobin protein (1-190AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A