

# CYP2E1 Antibody, HRP conjugated



PACO26102

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Metabolizes several precarcinogens, drugs, and solvents to reactive metabolites. Inactivates a number of drugs and xenobiotics and also bioactivates many xenobiotic substrates to their hepatotoxic or carcinogenic forms.

**Gene ID:**

CYP2E1

**Uniprot**

P05181

**Synonyms:**

Cytochrome P450 2E1 (EC 1.14.13. -) (4-nitrophenol 2-hydroxylase) (EC 1.14.13. n7) (CYP11E1) (Cytochrome P450-J), CYP2E1, CYP2E

**Immunogen:**

Recombinant Human Cytochrome P450 2E1 protein (49-288AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A